## LINET

# MULTICARE X ACCESSORIES





























#### **RECOMMENDED KIT**

FOLDABLE IV & DRIVE

VENTILATOR CIRCUIT HOLDER

KIT OF STABILIZING PADS ADVANCED







#### LIFTING POLES & ACCESSORIES

The lifting pole is an accessory used to increase patient independence and mobility.



DESCRIPTION	SWL	SURFACE	MODEL
<b>Lifting pole</b> to be used with fixed head-end	75 kg	powder coated RAL 9002	11011410B0000
Lifting pole to be used with fixed head-end	75 kg	powder coated RAL 9006	11011410A0000
Lifting pole to be used with fixed head-end	75 kg	chrome plated	1101141000000



2

DESCRIPTION	SWL	SURFACE	MODEL
Grey plastic triangular holder to be placed on a lifting pole	75 kg	plastic grey color	4ROTGERSG700-2

LINET | Multicare X Accessories



DESCRIPTION	SWL	SURFACE	MODEL
Grey plastic triangular holder to be placed on a lifting pole adjustable by retractor	75 kg	plastic grey color	4ROTGNR100GR-2



DESCRIPTION	SWL	SURFACE	MODEL
IV holder, 3 hooks to be placed on a lifting		powder coated RAL 9002	16010700B0000
pole		powder coated RAL 9006	16010700A0000



DESCRIPTION	SWL	SURFACE	MODEL
IV holder, 3 hooks to be placed on a lifting pole		stainless steel	4MAPL00N1001

# IV POLES & ACCESSORIES

Telescopic IV poles can be inserted into the sockets or attached on the bed frame at the head end.



DESCRIPTION	SWL	SURFACE	MODEL
Foldable IV & Drive – telescopic infusion stands and driving handles with i-Drive Power controller holder. Intended for use without head board. (Attachable on bed frame under head end.)	10 kg (5 kg per hook)	powder coated	11022200Y0MCX

SWL = Safe Working Load



DESCRIPTION	SWL	SURFACE	MODEL
Foldable IV & Drive – telescopic infusion stands and driving handles. Intended for use without head board. (Attachable on bed frame under head end.)	10 kg (5 kg per hook)	powder coated	11022200N0MCX



DESCRIPTION	SWL	SURFACE	MODEL
IV pole, with plastic box, 4 hooks, height adjustable, bush size D19	8 kg (2 kg per hook)	chrome plated	4PKV26107200



DESCRIPTION	SWL	SURFACE	MODEL
Infusion bottle basket to be hang on IV pole	Defined by SWL of IV pole	stainless steel	4DR426101

## LATERAL TILT SUPPORT

Recommended accessories for making lateralization even safer.



DESCRIPTION	SWL	SURFACE	MODEL
Ventilator circuit holder (if ordered separatelly - not as the part of KIT 3)	0.6 kg	stainless steel, plastic	4MAS25012600

LINET | Multicare X Accessories



DESCRIPTION	SWL	SURFACE	MODEL
Kit of stabilizing pads ADVANCED for the lateral tilt (6pcs foam frame pads + 15 pcs pearl balls fixation pillows)	_	mattress foam and pearls / washable cover	4ZE104041317



DESCRIPTION	SWL	SURFACE	MODEL
Kit of stabilizing pads for the lateral tilt (if ordered separatelly - not as the part of KIT 3)	_	mattress foam / washable	4GMKLS000003A



DESCRIPTION	SWL	SURFACE	MODEL
Positioning cushion for lateral tilt, washable stretch cover (1 piece)	_	mattress foam / washable	4SET17500001

### HIGHER SAFETY & COMFORT

Enhance further safety and comfort of the patient.



DESCRIPTION	SWL	SURFACE	MODEL
Pair of Protectors, complete set with cover for Multicare	-	powder coated RAL 9006 / HPL	11005900B0MCR-1

SWL = Safe Working Load 5



DESCRIPTION	SWL	SURFACE	MODEL
m-Panel - displays statuses of the bed (compatible only with MPA:01)	-		4\$6018690



DESCRIPTION	SWL	SURFACE	MODEL
Handset holder for Multicare	_	plastic	1102950000000

### **SHELVES**

Shelves and tray to be placed at the bed end for further medical equipment.



DESCRIPTION	SWL	SURFACE	MODEL
<b>Monitor shelf</b> for Multicare, desk dimensions 620 × 330 mm	15 kg	plastic / powder coated RAL 9006	11026300A0006



6

DESCRIPTION	SWL	SURFACE	MODEL
Writting shelf for plastic or metal bed ends	10 kg	plastic / powder coated RAL 9002	11023230B0000

LINET | Multicare X Accessories



DESCRIPTION	SWL	SURFACE	MODEL
Universal utility shelf	10 kg	powder coated RAL 9002	11025800B0000

### **OXYGEN BOTTLE HOLDERS**

Oxygen bottle holders can also be fitted to the bed frame. Multicare X provides 3 separate oxygen bottle mounting positions.



DESCRIPTION	SWL	SURFACE	MODEL
Horizontal oxygen bottle (5I) holder to be placed on crossbar on the head end of the bed (Size of O <sub>2</sub> bottle Max: L520 × D140; Min. L290 × D90)	15 kg	stainless steel	4MAR2010PC004



DESCRIPTION	SWL	SURFACE	MODEL
Vertical oxygen bottle holder to be placed on head end (Size of $O_2$ bottle: Max: L800 $\times$ D140; Min: L325 $\times$ D120)	15 kg	stainless steel	4MAR2010PC007-1



DESCRIPTION	SWL	SURFACE	MODEL
Vertical oxygen bottle holder to be placed on head end (Size of O <sub>2</sub> bottle Max: L800 × D120; Min: L365 × D100)	15 kg	stainless steel	4MAR2015PC001

SWL = Safe Working Load 7



DESCRIPTION	SWL	SURFACE	MODEL
Adaptor for the O <sub>2</sub> bottle holder 4MAR2010PC007-1	25 kg	powder coated RAL 9002	11029700A0000

### URINE/DRAINAGE BOTTLES & SACK HOLDERS

The stainless steel rail and undercarriage provide hooks and sockets that can be used for urine / drainage bags and bottles, attaching cannulas and more.



DESCRIPTION	SWL	SURFACE	MODEL
Urine bottle / shoes holder with two hooks	5 kg	powder coated RAL 9002	11023800B0000



DESCRIPTION	SWL	SURFACE	MODEL
Urine sack holder, level self-adjustable	5 kg	powder coated RAL 9002	11023700B0000
	5 kg	powder coated RAL 9006	11023700A0000



DESCRIPTION	SWL	SURFACE	MODEL
Urine sack holder, on undecarriage	_	powder coated RAL 9002	11023820B0000





